

Tarse

Upmarket podiatry chair with very stable structure and high comfort. Specially designed to optimize workspace due to chair 240° rotation. 3 motors control height, backrest and seat inclination. Independently adjustable legrests by means of gas lift. Reversible arms and adjustable headrest. Also available 5 motors model (2246C) to control both legrests movement. White PU upholstery of high quality and easy cleaning.

Electric Podiatry Chair

Ref. 2246A

Size without armrest: 186x57 cm

Size with armrest: 186x82 cm

Weight: 79 Kg

Upholstery: PU

Motors: 3

1 Breather hole

2 Angle

3 Extension ±8 cm

4 Angle 0° ~ 80°

5 Angle 0° ~ 10°

6 Rotation 0° ~ 240°

7 Angle 0° ~ 85°

8 Up and Down 56/89 cm

9 Extension ±14 cm

Electric Podiatry Chair

Ref. 2246C

Size without armrest: 186x57 cm

Size with armrest: 186x82 cm

Weight: 79 Kg

Upholstery: PU

Motors: 5

OPTIONS

Bed cover

10 Foot controller Ref. FC-003/FC-005

Wheels

Heating

Wireless

Memory

Silicone

11 Debris tray Ref. 30030

PACKING SIZE cm

136x84x64

N.W. kg

79

G.W. kg




92

MAX WEIGHT

180 kg



Tarse



Tarse

